

Electrostatic methods:
EFM, SPFM, KPFM (KFM)

The Kelvin method

In this method we measure the current induced in the wire connecting two materials of difference work functions when the distance between the two is varied in a periodic way. The current is due to the change in capacitance and the contact potential difference:

$$Q = C.V = C(z).(\phi_1 - \phi_2)$$

[$= \epsilon S/z.(\phi_1 - \phi_2)$ for a parallel plate capacitor]

Modulating the distance z , for example with a sine wave: $z(t) = z_0 + z_1 \sin(\omega t)$

$$I(t) = \partial Q / \partial t = \partial C / \partial z. z_1 \omega \cos(\omega t).(\phi_1 - \phi_2) = I_0 \cos(\omega t)$$



The measurement consists in applying a voltage difference between the two such that the current I is zero:

$$I_0 = 0 = \partial C / \partial z. z_1.(\phi_1 - \phi_2 + \Delta V)$$



That gives $\phi_1 - \phi_2 = -\Delta V$

We will later see how this technique can be implemented in an Atomic Force Microscope. In this way we will be able to “map out” the local value of $\phi(x,y)$ with very high spatial resolution (nm).



Electrostatic forces: Scanning Polarization Force Microscopy



Non-contact image




Contact Image



c)

Dielectric Spectroscopy

Freq. Dependence of $\epsilon(\omega)$



Modulation techniques in SPFM

Electrostatic force for a dielectric material

$$F = -\frac{1}{2} \times \frac{\partial C}{\partial z} \times (V - \varphi)^2 = \tilde{\varepsilon}(\omega) \times (V - \varphi)^2$$

Modulating $V = V_0 + V_1 \sin \omega t$


$$F_{dc} = \tilde{\varepsilon}(\omega) \times f(z) \times (V_0 - \varphi)^2 \quad \text{Topography + surface potential}$$

$$F_{\omega} = \tilde{\varepsilon}(\omega) \times f(z) \times (V_0 - \varphi) \times V_1 \quad \text{Surface potential}$$

$$F_{2\omega} = \tilde{\varepsilon}(\omega) \times f(z) \times V_1^2 \quad \longrightarrow \quad \text{Topography}$$

Measuring contact potentials: Kelvin Probe Microscopy



Liquid crystal 8CB:




Topography

Surface potential




Original isotropic drop





Smectic bilayers on top of a monolayer




Kelvin Probe Force Microscopy



Kelvin Probe Force Microscopy



Surface Potential Imaging of PbS






Electron
accumulation

Hole
accumulation

topography

SP $V_g=52V$

SP $V_g=-52V$



0 nm

10.6 nm



4.80 V

4.75 V



5.07 V

5.02 V